



# Rapid Phenotypic Detection of Microbial Resistance in Gram-Positive Bacteria by a Real-Time Laser Scattering Method<sup>†</sup>

Evgeny A. Idelevich<sup>1</sup>, Matthias Hoy<sup>1</sup>, Dennis Görlich<sup>2</sup>, Dennis Knaack<sup>1</sup>, Barbara Grünastel<sup>1</sup>, Georg Peters<sup>1</sup>, Matthias Borowski<sup>2‡</sup> and Karsten Becker<sup>1</sup>\*

<sup>1</sup> Institute of Medical Microbiology, University Hospital Münster, Münster, Germany, <sup>2</sup> Institute of Biostatistics and Clinical Research, University of Münster, Münster, Germany

### **OPEN ACCESS**

### Edited by:

Miklos Fuzi, Semmelweis University, Hungary

### Reviewed by:

Sebastian Guenther, Freie Universität Berlin, Germany César de la Fuente, Massachusetts Institute of Technology, United States Osmar Nascimento Silva, Universidade Católica Dom Bosco, Brazil

### \*Correspondence:

Karsten Becker kbecker@uni-muenster.de <sup>‡</sup>These authors have contributed equally to this work.

### Specialty section:

This article was submitted to Antimicrobials, Resistance and Chemotherapy, a section of the journal Frontiers in Microbiology

**Received:** 23 March 2017 **Accepted:** 29 May 2017 **Published:** 14 June 2017

### Citation:

Idelevich EA, Hoy M, Görlich D, Knaack D, Grünastel B, Peters G, Borowski M and Becker K (2017) Rapid Phenotypic Detection of Microbial Resistance in Gram-Positive Bacteria by a Real-Time Laser Scattering Method. Front. Microbiol. 8:1064. doi: 10.3389/fmicb.2017.01064 We developed a methodology for antimicrobial susceptibility testing (AST) based on the BacterioScan<sup>TM</sup>216R laser scattering technology, using methicillin resistance in Staphylococcus aureus and vancomycin resistance in enterococci as exemplar for important resistance phenotypes. Fifty methicillin-resistant (MRSA) and 50 methicillin-susceptible (MSSA) S. aureus, as well as 50 vancomycin-resistant enterococci (VRE) and 50 vancomycin-susceptible enterococci (VSE) isolates were used for the study. Optimal test conditions were derived by investigating the effects of inoculum size, medium, incubation temperature and broth filtration. We proposed four different statistical approaches for rapid discrimination between resistant and susceptible bacteria. The statistical approach based on raw measurements of bacterial concentrations delivered sensitivity of 100% and specificity of 94% for discrimination between MRSA and MSSA already after 3 hours of incubation. Categorical agreement of ≥90% was achieved after 140 min with this approach. Differentiation between VRE and VSE was possible with 98% sensitivity and 92% specificity after 3 hours, using a sophisticated statistical approach based on concentration slopes derived from the raw concentration measurements. This approach provided categorical agreement of  $\geq$ 90% after 165 min. The sensitivity and specificity estimates were confirmed by leave-one-out cross validation. In conclusion, the phenotypic AST methods developed in this study are promising for rapid detection of MRSA and VRE. The development and application of this technology would allow early detection of the resistant pathogens, thus facilitating swift change to the targeted antimicrobial treatment as well as timely initiation of appropriate infection control measures. Further studies are warranted to validate this approach for the detection of other resistance phenotypes, including direct testing from clinical specimens.

Keywords: rapid diagnostics, resistance detection, susceptibility, laser scattering, BacterioScan

<sup>†</sup>This work was presented in part at the 26th European Congress of Clinical Microbiology and Infectious Diseases (ECCMID), 9th – 12th April 2016, Amsterdam, Netherlands [P0847]).

# INTRODUCTION

Increasing prevalence of multidrug resistant microorganisms (MDROs) emphasizes the need for their rapid detection and identification (Roca et al., 2015). Early knowledge of pathogen's species and antimicrobial susceptibility profile allows prompt change from empirically chosen antimicrobial treatment to the appropriate and targeted therapy (Doern et al., 1994; Idelevich et al., 2015). Administration of specific antibiotics not only enables effective treatment but also allows to minimize collateral damage on physiological flora and to prevent the development of resistant bacteria (Goldstein, 2011). Furthermore, rapid detection of MDROs facilitates swift initiation of infection control measures (Roca et al., 2015). While virtually immediate identification of microbial cultures has become available after the introduction of matrix-assisted laser desorption ionization time-of-flight mass spectrometry (MALDI-TOF MS) in clinical microbiology (Bizzini and Greub, 2010), antimicrobial susceptibility testing (AST) still takes long time (van Belkum et al., 2013). In diagnostic routine, the results of AST are usually available only on the next day or later after test initiation even though optimized protocols may accelerate AST (Kerremans et al., 2008; Idelevich et al., 2014a,b). Molecular methods can provide information on the presence of resistance genes more rapidly, but they are expensive and do not always correspond with the phenotypic testing (Jorgensen and Ferraro, 2009). Thus, rapid methods for detection of resistant phenotypes are urgently needed.

In this study, we focused on methicillin resistance in *S. aureus* and vancomycin resistance in enterococci as exemplar for resistance phenotypes with enormous medical and socio-economic burden (Köck et al., 2010). Methicillin-resistant *Staphylococcus aureus* (MRSA) is a major nosocomial pathogen and its rapid differentiation from methicillin-susceptible *S. aureus* (MSSA) isolates is crucial for the initiation of adequate infection control measures and antimicrobial therapy (Köck et al., 2014). Vancomycin-resistant enterococci (VRE) represent another multidrug resistant organism with increasing prevalence and similar challenges for the health care system (Gastmeier et al., 2014; Humphreys, 2014).

Here, we aimed at developing methods for AST that are based on the laser scattering technology implemented in the BacterioScan<sup>TM</sup>216R (BacterioScan Inc., St Louis, MO, United States) instrument. BacterioScan uses the laser scattering technology to quantify bacteria in up to 16 liquid samples in real time. The combination of the laser scattering technology with sophisticated statistical real-time analyses may provide a promising approach for fast and reliable discrimination between susceptible and resistant strains. We propose four statistical approaches for this purpose. The study's main objectives were (i) to quantify the discriminatory performance of each method in terms of sensitivity, specificity, and categorical agreement and (ii) to determine the incubation time required for a reliable discrimination between resistant and susceptible isolates.

# MATERIALS AND METHODS

# **Investigation of Optimal Test Conditions**

Comparative investigation of optimal conditions for AST by laser scattering included the evaluation of the effects of inoculum size, various media, broth filtration as well as the effect of incubation temperature. All experiments were performed in triplicate.

### Effect of Inoculum Size

After preparation of bacterial suspensions of reference strains S. aureus ATCC BAA-44 and Escherichia coli ATCC 25922 in brain-heart infusion (BHI) broth (Merck, Darmstadt, Germany) with 0.5 McFarland turbidity using a nephelometer (Densimat, bioMérieux, Marcy l'Etoile, France), serial dilutions in BHI were made to produce suspensions with estimated 10-fold concentrations from 1 cfu/ml to  $1 \times 10^8$  cfu/ml. Additionally, a suspension with  $5 \times 10^5$  cfu/ml was prepared, which reflects the bacterial density recommended as starting inoculum for AST by the International Organization for Standardization ISO (2006) and the Clinical Laboratory Standards Institute CLSI (2015) guidelines. The real cell concentration in these suspensions was confirmed by vital cell counting after plating of a sample onto tryptic soy agar (TSA) plates in triplicate. Two ml samples of each dilution and a sterile control sample were added to the cuvettes of the BacterioScan instrument, and the bacterial concentration was measured for 24 h at 35°C. Applying different starting inocula, the time until fivefold and 10-fold increase in inoculum was documented to determine the length of time needed to detect the growth by laser scattering.

### Effect of Medium

Bacterial suspensions with 0.5 McFarland turbidity (Densimat, bioMérieux) were prepared in cation-adjusted Mueller-Hinton broth (CA-MHB, BD Diagnostics, Heidelberg, Germany), BHI broth (Merck), tryptic soy broth (TSB, BD Diagnostics, Germany) and lysogeny broth (LB, BD Diagnostics, Germany) for *S. aureus* ATCC BAA-44 and *Escherichia coli* ATCC 25922. All broths were used unfiltered. The suspensions were diluted 1:200 to produce an inoculum concentration of approximately  $5 \times 10^5$  cfu/ml. The real cell concentration was verified by vital cell count after plating onto TSA plates. Two ml samples were measured in BacterioScan for 24 h at 35°C. The time until fivefold and 10-fold concentration increase was recorded for comparison of the length of times needed to detect the growth using different media.

### Effect of Broth Filtration

Suspensions of 0.5 McFarland turbidity (Densimat, bioMérieux) were prepared in CA-MHB and BHI broth for *S. aureus* ATCC BAA-44 and *Escherichia coli* ATCC 25922. One sample set was filtered through a 0.2  $\mu$ l syringe filter (Sarstedt, Nümbrecht, Germany) while other sample set was left unfiltered. The suspensions were diluted 1:200 to produce an inoculum concentration of approximately 5  $\times$  10<sup>5</sup> cfu/ml (confirmed by vital cell count). Two ml samples were incubated in BacterioScan for 8 h at 35°C. The time until fivefold and 10-fold

increase in inoculum was documented to compare the length of time needed for growth detection using filtered or unfiltered broth.

### Effect of Incubation Temperature

Applying the BacterioScan instrument, the effect of incubation temperature on staphylococcal AST was investigated on MSSA reference strains ATCC 29213 and ATCC 25923 as well as MRSA reference strains ATCC BAA-44 and ATCC 43300 using CA-MHB as medium. One ml of cefoxitin solution was added to 1 ml of bacterial suspension to produce a final inoculum of  $5 \times 10^5$  cfu/ml and a final cefoxitin concentration of 4 µg/ml (breakpoint concentration of cefoxitin, which allows the differentiation between MRSA and MSSA, according to the European Committee on Antimicrobial Susceptibility Testing EUCAST (2016) and CLSI (2016). The temperatures investigated were 30°C, 33°C, 34°C, 35°C, 36°C, and 37°C (Supplementary Figure S1). Real microbial concentrations were confirmed by vital cell count. Controls without antibiotics were implemented.

## **Bacterial Strains**

Fifty MRSA and 50 MSSA isolates, as well as 50 vancomycin-resistant *E. faecium* (VR-*E. faecium*) and 50 vancomycin-susceptible *E. faecium* (VS-*E. faecium*) clinical isolates were used in the study. Only one isolate per patient was included. Species identification was confirmed by MALDI-TOF MS. Methicillin resistance in *S. aureus* and vancomycin resistance in *E. faecium* isolates were confirmed by GenoType MRSA (Hain Lifescience, Nehren, Germany) assay and *vanA/vanB* PCR, respectively. Among VREs, 24 and 26 isolates were *vanA*-positive and *vanB*-positive, respectively.

### Antimicrobial Susceptibility Testing Reference Method

Minimum inhibitory concentrations (MIC) of cefoxitin for *S. aureus* and vancomycin for *E. faecium* were determined by broth microdilution (BMD) reference method according to the ISO 20776-1 (2006) and CLSI (2015) guidelines. Cefoxitin and vancomycin powders were purchased from Sigma-Aldrich (Saint Louis, MO, United States). The range of tested concentrations was 0.25–128  $\mu$ g/ml for both antimicrobials. The tests were performed in triplicate on different days and median MIC values were used for analysis. Reference strains *S. aureus* ATCC 29213 and *E. faecalis* ATCC 29212 were used on every testing day for staphylococci and enterococci, respectively.

# Measuring Microbial Concentration by Forward Laser Scattering

Quantitative estimation of the microbial concentration was performed using the BacterioScan<sup>TM</sup>216R instrument, which uses narrow angle forward laser scattering as a sensitive optical method for measuring light scattered by bacteria suspended in a liquid sample. Multiple simultaneous measurements of the forward scattering and optical density performed over time, in combination with proprietary calculation algorithms, provide accurate estimation of microbial concentrations even at low cfu/ml levels (Anbazhagan and Regelman, 2015). The tabletop instrument accommodates up to four disposable multi-cuvettes, each of them containing four wells of approximately 2 ml volume. Additionally to the possibility of incubation at room temperature, onboard incubation allows temperatures between  $30^{\circ}$ C und  $42^{\circ}$ C.

0.5 McFarland turbidity (Densimat) suspensions were prepared in CA-MHB (BD Diagnostics) and diluted to produce a starting inoculum concentration of approximately  $5 \times 10^5$  cfu/ml after adding the antibiotic. The cell concentration in these suspensions was confirmed by vital cell count after plating onto TSA plates in triplicate. The samples were incubated in the BacterioScan<sup>TM</sup>216R device for 6 h. For detection of methicillin resistance in staphylococci, cefoxitin was used in the breakpoint concentration of 4 µg/ml as recommended by CLSI (2016) and EUCAST (2016) and the samples were incubated at 34°C. For detection of vancomycin resistance in enterococci, the vancomycin was added to produce the final breakpoint concentration 4 µg/ml (CLSI, 2016; EUCAST, 2016), followed by the incubation of samples at 35°C. Antibioticfree growth controls were included for each isolate. The measurements were taken automatically approximately every 3 min for each sample, and then converted into sequences of one measurement per minute by simple linear interpolation. Example growth curves are given in **Figure 1**. While the growth of susceptible strains is inhibited by antibiotic, the growth curves with a given antibiotic compound are comparable to those of the growth control without antibiotic in resistant strains.

# **Statistical Analysis**

### **Optimal Test Conditions**

The effects of inoculum size, various media, and broth filtration were analyzed by univariate ANOVA with Tukey's test using IBM SPSS Statistics 23 software (IBM, Armonk, NY, United States).

# Approaches to Discriminate between Resistant and Susceptible Bacteria

Four different approaches to discriminate between resistant and susceptible bacteria based on data delivered by the laser scattering methodology were investigated (Supplementary Figure S2) based on

- (a) concentration,
- (b) ratio of concentrations (sample with antibiotics vs. growth control),
- (c) concentration slope, and
- (d) ratio of concentration slopes (sample with antibiotics vs. growth control).

In approach (a) the concentration measurements  $x_t$  at time t, say t = 120 [minutes], of the 50 resistant and 50 susceptible isolates (each with added antibiotic) were used to perform a receiver-operator-curve (ROC) analysis, in order to determine the discrimination cutoff at time t = 120 that delivered the highest sum of sensitivity and specificity. This was done for all times  $t = 60, \dots, 360$ . Approach (b) was equal to (a), except that the concentration ratio  $x_t/y_t$  instead of the concentration



 $x_t$  was considered for ROC analyses, where  $y_t$  denotes the concentration of the growth control that corresponds to  $x_t$  In approach (c), the slope of bacterial growth  $s_t$  was derived from the time series of concentration measurements  $x_t$  at time t = 60, ..., 360 by means of the SCARM filter (Borowski and Fried, 2014). The ROC analyses were then performed on  $s_t$  at time t = 60, ..., 360. The SCARM is a real-time filter, i.e., it provides a slope estimate  $s_t$  as soon as the concentration  $x_t$  has been measured. Approach (c) could therefore be applied in practice. The rationale for choosing this approach as well as more details about the SCARM estimation of  $s_t$  are given in the Supplementary Description. In approach (d), the ratio  $s_t/u_t$  was considered for ROC analysis, where  $u_t$  denotes the SCARM slope estimate of the growth control  $Y_t$  that corresponds to  $x_t$ 

For each approach (a)-(d), the ROC analyses delivered a best discrimination cutoff and its sensitivity and specificity estimates (with 95% confidence intervals) at time t = 60, ..., 360. Additionally, the incubation duration until overall categorical agreement (CA) reached  $\geq$ 90%, which is an acceptance limit for accuracy of susceptibility testing recommended by ISO (2007) 20776-2 guideline, was calculated. Since cutoff determination and estimation of sensitivity and specificity were carried out on the same data, probably leading to biased estimates, we

calculated all sensitivity and specificity estimates and their 95% confidence intervals by leave-one-out cross-validation, too. All analyses described here were performed using R, version 3.3.2.

# RESULTS

# **Investigation of Optimal Test Conditions**

For both reference strains, no significant difference between starting inocula from  $5 \times 10^5$  to  $1 \times 10^8$  was demonstrated regarding time to detection of fivefold or 10-fold increase in bacterial concentration. However, this time was significantly longer for lower starting inocula (Supplementary Table S1). Neither the four liquid media nor the broth filtration influenced the time to growth detection significantly (Supplementary Tables S2, S3). The MRSA reference strain ATCC BAA-44 grew, as expected, in the presence of cefoxitin at incubation temperatures 30°C, 33°C, 34°C, 35°C, 36°C. However, at 37°C it was completely inhibited by cefoxitin, thus behaving like an MSSA (Supplementary Figure S1). The growth of the other MRSA reference strain ATCC 43300 was also delayed to some extent at 37°C, but the growth was restored after approximately 12 hours (data not shown). The MSSA reference strains ATCC 29213 and ATCC 25923 were inhibited by cefoxitin at all incubation temperatures, as expected (data not shown).

## **Standard Resistance Detection**

MIC<sub>50</sub>, MIC<sub>90</sub> and MIC range, as determined by broth microdilution, were 4  $\mu$ g/ml, 4  $\mu$ g/ml, 2–4  $\mu$ g/ml, and 64  $\mu$ g/ml,  $\geq$ 256  $\mu$ g/ml, 16– $\geq$ 256  $\mu$ g/ml for MSSA and MRSA, respectively. For VS-*E. faecium* and VR-*E. faecium*, MIC<sub>50</sub>, MIC<sub>90</sub> and MIC ranges were 1  $\mu$ g/ml, 1  $\mu$ g/ml, 0.5 – 2  $\mu$ g/ml and  $\geq$ 256  $\mu$ g/ml,  $\geq$ 256  $\mu$ g/ml, 32– $\geq$ 256  $\mu$ g/ml, respectively.

# Performance of the Approaches to Discriminate between Resistant and Susceptible Bacteria

Estimates of the sensitivity and specificity of the approaches (a)–(d) applied to discriminate between MRSA and MSSA are presented in **Table 1**. For all approaches, both the sensitivity and specificity were >90% after 4 h, and >98% after 5 h. Interestingly, the approaches (a) and (c), i.e., those approaches that do not require growth controls, performed better than approaches (b) and (d). For approach (a) and (c), the estimates of sensitivity and specificity were >90% already after 3 h, with lower bounds of the 95% CIs of >89.4% for sensitivity and >83.5% for specificity. The overall categorical agreement of  $\geq$ 90% was achieved after 140 and 147 min with approaches (a) and (c), respectively, followed by 151 and 166 min with approaches (b) and (d), respectively.

**Table 2** shows the sensitivity and specificity estimates for VR-*E. faecium* and VS-*E. faecium*. The approaches (c) and (d), i.e., the approaches based on concentration slopes, performed best, with slightly larger values for approach (c). This approach delivered sensitivity and specificity of >90% after 3 h, where the lower bounds of the 95% CIs were >89.4% for sensitivity and >80.8 for specificity. While approaches (c) and (d) provided categorical agreement  $\geq$ 90% after 165 min and 164 min, respectively, this limit was not achieved with approach (a) before 215 min, and was not achieved at all with approach (b).

The sensitivity and specificity estimates and their 95% CIs that were obtained by leave-one-out cross-validation (Supplementary Tables S4, S5) were very close to those presented in **Tables 1**, **2**, indicating their reliability.

# DISCUSSION

Rapid AST is a function of several components addressed in our study: optimal testing conditions for rapid and correct exhibiting of phenotypic resistance, susceptible growth detection technology, and sophisticated statistical algorithms for early differentiation between resistant and susceptible phenotypes.

Our investigations of the effects of different testing conditions for the optimal, i.e., rapid and exact, AST revealed that it was unnecessary to change conditions which are recommended by ISO (2006) and CLSI (2015) guidelines for reference broth microdilution method. Use of other broths than recommended CA-MHB did not accelerate the detection of

Approach		Ŭ	a)			U	(q			J	(1)			(q)		
Min.	SE (%)	(95% CI)	SP (%)	(95% CI)	SE (%)	(95% CI)	SP (%)	(95% CI)	SE (%)	(95% CI)	SP (%)	(95% CI)	SE (%)	(95% CI)	SP (%)	(95% CI)
60	84	(70.9–92.8)	48	(33.7-62.6)	52	(37.4–66.3)	78	(64-88.5)	66	(51.2–78.8)	64	(49.2–77.1)	48	(33.7-62.6)	76	(61.8–86.9)
120	82	(68.6–91.4)	80	(66.3–90)	80	(66.3–90)	88	(75.7–95.5)	78	(64 - 88.5)	84	(70.9–92.8)	68	(53.3-80.5)	06	(78.2–96.7)
180	100	(92.9–100)	94	(83.5–98.7)	86	(73.3–94.2)	98	(89.4–99.9)	98	(89.4–99.9)	96	(86.3–99.5)	84	(70.9–92.8)	100	(92.9-100)
240	100	(92.9–100)	100	(92.9–100)	94	(83.5-98.7)	94	(83.5–98.7)	100	(92.9–100)	100	(92.9–100)	100	(92.9–100)	100	(92.9-100)
300	100	(92.9–100)	100	(92.9–100)	100	(92.9–100)	100	(92.9–100)	100	(92.9–100)	100	(92.9–100)	100	(92.9–100)	100	(92.9-100)
360	100	(91.6-100)	98	(89.4–99.9)	100	(91.6-100)	98	(89.4–99.9)	100	(91.6–100)	98	(89.4–99.9)	100	(91.6–100)	100	(92.9-100)

Approach		-	a)			~	(q			J	(;			(q)		
Min.	SE (%)	(95% CI)	SP (%)	(95% CI)	SE (%)	(95% CI)	SP (%)	(95% CI)	SE (%)	(95% CI)	SP (%)	(95% CI)	SE (%)	(95% CI)	SP (%)	(95% CI)
00	06	(78.2–96.7)	30	(17.9-44.6)	06	(78.2-96.7)	44	(30-58.7)	56	(41.3–70)	62	(47.2–75.3)	48	(33.7-62.6)	78	(64-88.5)
120	82	(68.6–91.4)	66	(51.2–78.8)	78	(64-88.5)	84	(70.9–92.8)	86	(73.3–94.2)	70	(55.4–82.1)	80	(66.3–90)	82	(68.6–91.4)
180	92	(80.8–97.8)	86	(73.3–94.2)	06	(78.2–96.7)	82	(68.6–91.4)	98	(89.4–99.9)	92	(80.8–97.8)	94	(83.5–98.7)	98	(89.4–99.9)
240	92	(80.8–97.8)	06	(78.2–96.7)	92	(80.8–97.8)	80	(66.3–90)	98	(89.4–99.9)	92	(80.8–97.8)	98	(89.4–99.9)	94	(83.5–98.7)
300	06	(78.2–96.7)	98	(89.4–99.9)	96	(86.3–99.5)	76	(61.8–86.9)	100	(92.9–100)	98	(89.4–99.9)	100	(92.9–100)	92	(80.8–97.8)
360	06	(78.2–96.7)	100	(92.9–100)	96	(86.3-99.5)	78	(64–88.5)	98	(89.4–99.9)	96	(86.3–99.5)	98	(89.4–99.9)	86	(73.3–94.2)

Idelevich et al

(Supplementary Table S1). Broth filtration might be deemed helpful for testing by optical systems, but it also did not influence the time to growth detection considerably. Thus, we demonstrated that there is generally no reason to alter the BacterioScan testing conditions compared to those recommended for the reference method. On the other hand, it should be kept in mind that the introduction of deviations in test conditions might affect the test results. This can be best exemplified by the effect of different incubation temperatures. Elevation of incubation temperature up to only 37°C had a tremendous effect on the results of staphylococcal AST in our experiments. The MRSA strain ATCC BAA-44, which grew as expected in the presence of breakpoint cefoxitin concentration at temperatures from 30 to 36°C, was inhibited with the same antibiotic concentration at 37°C, i.e., it behaved as MSSA at higher temperature (Supplementary Figure S1). This effect was observed in MRSA strain ATCC BAA-44, but not in MRSA strain 43300, demonstrating the strain-dependence of this effect. The loss of methicillin resistance in some S. aureus strains at higher temperature is known and was described early by other investigators (Asheshov, 1966). Interestingly, such findings provide experimental background for the hypothesis that β-lactam antibiotics can be used in treatment of infections caused

with non- $\beta$ -lactam antibiotics. Clinical studies investigating the practical benefit of such regimens are warranted (Dhand et al., 2011; Poulsen et al., 2013). Another example of unfavorable effects of deviating test conditions from the standard is the so-called inoculum effect well described for some groups of antimicrobials, particularly for  $\beta$ -lactam antibiotics (Sabath et al., 1975; Brook, 1989). Therefore, unnecessary deviations in test conditions should be strictly avoided. As a result of the preliminary experiments and the theoretical considerations, we decided to keep the test conditions for the main experiments (broth type, inoculum size, antibiotic concentration) as close as possible to the official

by MRSA, e.g., in the framework of combination regimens

size, antibiotic concentration) as close as possible to the official recommendations for AST (ISO, 2006; CLSI, 2016; EUCAST, 2016). Only the sample volume (2 ml) had to be adapted to the requirements of the test device.

The key difference and concurrent major advantage of the approaches investigated in this study was the earlier evaluation of growth as compared to the standard testing methods, i.e., the visual reading of turbidity after 16–20 h with the reference broth microdilution method (ISO, 2006; CLSI, 2016). This was accomplished using the technology of cell counting implemented in the BacterioScan instrument in combination with the statistical algorithms for early differentiation between multiplying or inhibited microorganisms.

Reliable differentiation between MSSA and MRSA and between VSE and VRE was possible already after less than

growth (Supplementary Table S2). Similarly, inoculum sizes higher than recommended  $5 \times 10^5$  cfu/ml did not considerably accelerate the detection of growth. Merely, the inocula lower than  $5 \times 10^5$  cfu/ml resulted in the significantly delayed growth detection, as expected, because they first need to reach the instrument's lower detection limit during the incubation 3 h (**Tables 1**, **2**). Results from leave-one-out cross-validation were similar, indicating that the proposed methods would be as accurate and fast when applied to further concentration measurements of staphylococci and enterococci.

Attempts were made to find out the optimal method of statistical analysis, which allows for discriminating between susceptible and resistant isolates as early as possible. Interestingly, incorporating the growth control data into the statistical approaches [approaches (b) and (d)] resulted in neither faster nor more accurate differentiation between susceptible and resistant isolates of both staphylococci and enterococci. It even reduced accuracy in some situations, which might be due to the additional variability coming from each isolate's growth control. While testing staphylococci, the slope estimating approaches (c) and (d) had no advantage over the simpler approaches (a) and (b) that are based on concentrations only. However, the slope estimating approaches provided more accurate and faster discrimination between VRE and VSE. This is most probably because the concentrations of MRSA and MSSA differed considerably after 3 to 4 hours, whereas the concentrations of VRE and VSE did not (Supplementary Description). However, the slopes of VRE and VSE did differ considerably quite early, therefore making the slope estimating approaches (c) and (d) superior to the simpler "concentration approaches" (a) and (b) for enterococci. Since such structures are expectable also in concentration time series of other bacteria, we think that the slope estimating approaches can be more promising for some microorganism groups. Moreover, we are convinced that refining and combining statistical approaches could improve growth-based AST even further.

Early studies have identified light scattering technology as a promising tool for detecting microbial resistance within a short time (Berkman et al., 1970; Wyatt, 1972; Murray et al., 1980). Already at that time, efforts were made to develop instruments for rapid AST (Phillips, 1971; Stull, 1973). BacterioScan device has been recently introduced as an instrument for microbial quantification using laser scattering (Anbazhagan and Regelman, 2015). The device has an advantage of compact design and easy-to-use cuvettes, which make it promising particularly for the on-demand single resistance tests. Two very recent studies have investigated laser scattering technology using BacterioScan instrument for AST (Bugrysheva et al., 2016; Hayden et al., 2016). Bugrysheva et al. have demonstrated that laser scattering technology reduced the time to AST result by 50 to 75% for biothreat bacteria Bacillus anthracis, Yersinia pestis, and Burkholderia pseudomallei, as compared to conventional methods (Bugrysheva et al., 2016). Bugrysheva et al. (2016) discriminate between growing resistant and susceptible bacteria by comparing CFU/ml counts for a bacterial strain with and without an antimicrobial within the 30-min intervals by means of *p*-values from paired *t*-tests and Wilcoxon-tests. Hayden et al. (2016) have shown accurate determination of antimicrobial susceptibility using three isolates each of S. aureus, E. coli, and Pseudomonas aeruginosa. Hayden et al. (2016) consider rules based on areas under the curve of laser scattering optical density measurements to determine antimicrobial susceptibility. As a proof-of-principle, the authors

demonstrate a high degree of categorical concordance with two commercial methods for AST (Hayden et al., 2016). However, the quantification of reliability to discriminate between resistant and susceptible bacteria with the approaches used in these both studies is limited due to a low number of isolates used.

In our study, we directly used the bacterial concentrations calculated by the instrument for the analysis. We did not calibrate the calculated concentrations to the bacterial concentrations determined by the plating of suspensions for colony count, as it was done in another study (Bugrysheva et al., 2016). The statistical analysis was based on the concentration data directly generated by the instrument (calculated cfu/ml). The comparability of these primary data between resistant and susceptible isolates as well as samples with and without antibiotics is obvious because the proprietary software generates all data sets in the same way. In contrast, transforming those calculated concentrations to "calibrated cfu/ml" has no additional benefit for analysis and, to our opinion, can even result in additional bias. This is because the real concentrations determined by colony count method for a chosen "calibrating strain" does not necessarily reflect the concentrations of other isolates even within the same species. Furthermore, growth in chains, as pointed out by Bugrysheva et al. (2016), or in clusters of cohering cells may result in colony counting error. This could also apply to enterococci and staphylococci used in our study.

Since this study investigated the detection of single resistance phenotypes by using only one antibiotic against particular pathogen, the information on susceptibility toward other antibiotics is missing. Simultaneous testing of multiple antimicrobials would provide information on alternative treatment options and, thus, facilitate prompt and correct treatment choice.

# CONCLUSION

The phenotypic AST methods developed in this study are promising for rapid detection of MRSA and VRE. Further studies are required to investigate and validate these approaches for detecting other important resistance phenotypes. Future studies should also focus on direct testing from clinical specimens.

# **AUTHOR CONTRIBUTIONS**

EI and KB conceived and designed the experiments; MH, BG, and EI performed the experiments; EI, MH, DG, DK, GP, MB, and KB analyzed the data; MB and DG developed statistical algorithms and performed statistical analyses; EI wrote the manuscript, with input from MH, DG, DK, GP, MB, and KB.

# FUNDING

The BacterioScan instrument and consumables were provided free of charge for the study period by BacterioScan Inc.

# ACKNOWLEDGMENT

We are grateful to Damayanti Kaiser for expert technical assistance and to Rajesh Anbazhagan (BacterioScan Inc.) for helpful technical advices.

# REFERENCES

- Anbazhagan, R., and Regelman, D. V. (2015). "Novel real time bacterial growth and size measurement with laser scattering for improved antibiotic susceptibility testing," in *Poster 10/26. Abstracts of the 115th General Meeting of the ASM*, Washington, DC.
- Asheshov, E. H. (1966). Loss of antibiotic resistance in *Staphylococcus aureus* resulting from growth at high temperature. J. Gen. Microbiol. 42, 403–410. doi: 10.1099/00221287-42-3-403
- Berkman, R. M., Wyatt, P. J., and Phillips, D. T. (1970). Rapid detection of penicillin sensitivity in *Staphylococcus aureus*. *Nature* 228, 458–460. doi: 10.1038/228458a0
- Bizzini, A., and Greub, G. (2010). Matrix-assisted laser desorption ionization timeof-flight mass spectrometry, a revolution in clinical microbial identification. *Clin. Microbiol. Infect.* 16, 1614–1619. doi: 10.1111/j.1469-0691.2010. 03311.x
- Borowski, M., and Fried, R. (2014). Online signal extraction by robust regression in moving windows with data-adaptive width selection. *Statis. Comput.* 24, 597–613. doi: 10.1155/2011/143480
- Brook, I. (1989). Inoculum effect. *Rev. Infect Dis.* 11, 361–368. doi: 10.1093/clinids/ 11.3.361
- Bugrysheva, J. V., Lascols, C., Sue, D., and Weigel, L. M. (2016). Rapid antimicrobial susceptibility testing of *Bacillus anthracis*, *Yersinia pestis*, and *Burkholderia pseudomallei* by use of laser light scattering technology. *J. Clin. Microbiol.* 54, 1462–1471. doi: 10.1128/JCM. 03251-15
- CLSI (2015) Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; Approved Standard, 10th Edn. Wayne, PA: Clinical Laboratory Standards Institute.
- CLSI (2016) Performance Standards for Antimicrobial Susceptibility Testing, 26th Edn. Wayne, PA: Clinical Laboratory Standards Institute.
- Dhand, A., Bayer, A. S., Pogliano, J., Yang, S. J., Bolaris, M., Nizet, V., et al. (2011). Use of antistaphylococcal beta-lactams to increase daptomycin activity in eradicating persistent bacteremia due to methicillin-resistant *Staphylococcus aureus*: role of enhanced daptomycin binding. *Clin. Infect. Dis.* 53, 158–163. doi: 10.1093/cid/cir340
- Doern, G. V., Vautour, R., Gaudet, M., and Levy, B. (1994). Clinical impact of rapid in vitro susceptibility testing and bacterial identification. J. Clin. Microbiol. 32, 1757–1762.
- EUCAST (2016) Breakpoint Tables for Interpretation of MICs and Zone Diameters. The European Committee on Antimicrobial Susceptibility Testing. Version 6.0. Available at: www.eucast.org
- Gastmeier, P., Schröder, C., Behnke, M., Meyer, E., and Geffers, C. (2014). Dramatic increase in vancomycin-resistant enterococci in Germany. J. Antimicrob. Chemother. 69, 1660–1664. doi: 10.1093/jac/ dku035
- Goldstein, E. J. (2011). Beyond the target pathogen: ecological effects of the hospital formulary. *Curr. Opin. Infect. Dis.* 24(Suppl. 1), S21–S31. doi: 10.1097/01.qco. 0000393485.17894.4c
- Hayden, R. T., Clinton, L. K., Hewitt, C., Koyamatsu, T., Sun, Y., Jamison, G., et al. (2016). Rapid antimicrobial susceptibility testing using forward laser light scatter technology. *J. Clin. Microbiol.* 54, 2701–2706. doi: 10.1128/JCM.0 1475-16
- Humphreys, H. (2014). Controlling the spread of vancomycin-resistant enterococci. Is active screening worthwhile? J. Hosp. Infect. 88, 191–198. doi: 10.1016/j.jhin.2014.09.002
- Idelevich, E. A., Grunewald, C. M., Wüllenweber, J., and Becker, K. (2014a). Rapid identification and susceptibility testing of *Candida* spp. from positive blood cultures by combination of direct MALDI-TOF mass spectrometry and

# SUPPLEMENTARY MATERIAL

The Supplementary Material for this article can be found online at: http://journal.frontiersin.org/article/10.3389/fmicb. 2017.01064/full#supplementary-material

direct inoculation of Vitek 2. *PLoS ONE* 9:e114834. doi: 10.1371/journal.pone. 0114834

- Idelevich, E. A., Schüle, I., Grünastel, B., Wüllenweber, J., Peters, G., and Becker, K. (2014b). Acceleration of antimicrobial susceptibility testing of positive blood cultures by inoculation of Vitek 2 cards with briefly incubated solid medium cultures. J. Clin. Microbiol. 52, 4058–4062. doi: 10.1128/JCM. 02400-14
- Idelevich, E. A., Silling, G., Niederbracht, Y., Penner, H., Sauerland, M. C., Tafelski, S., et al. (2015). Impact of multiplex PCR on antimicrobial treatment in febrile neutropenia: a randomized controlled study. *Med. Microbiol. Immunol.* 204, 585–592. doi: 10.1007/s00430-014-0385-7
- ISO (2006). 20776-1. Clinical Laboratory Testing and In vitro Diagnostic Test Systems - Susceptibility Testing of Infectious Agents and Evaluation of Performance of Antimicrobial Susceptibility Test Devices - Part 1: Reference Method for Testing the In vitro Activity of Antimicrobial Agents Against Rapidly Growing Aerobic Bacteria Involved in Infectious Diseases. Geneva: International Organization for Standardization.
- ISO (2007). 20776-2. Clinical Laboratory Testing and In vitro Diagnostic Test Systems - Susceptibility Testing of Infectious Agents and Evaluation of Performance of Antimicrobial Susceptibility Test Devices - Part 2: Evaluation of Performance of Antimicrobial Susceptibility Test Devices. Geneva: International Organization for Standardization.
- Jorgensen, J. H., and Ferraro, M. J. (2009). Antimicrobial susceptibility testing: a review of general principles and contemporary practices. *Clin. Infect. Dis.* 49, 1749–1755. doi: 10.1086/647952
- Kerremans, J. J., Verboom, P., Stijnen, T., Hakkaart-van, R. L., Goessens, W., Verbrugh, H. A., et al. (2008). Rapid identification and antimicrobial susceptibility testing reduce antibiotic use and accelerate pathogen-directed antibiotic use. J. Antimicrob. Chemother. 61, 428–435. doi: 10.1093/jac/ dkm497
- Köck, R., Becker, K., Cookson, B., van Gemert-Pijnen, J. E., Harbarth, S., Kluytmans, J., et al. (2010). Methicillin-resistant *Staphylococcus aureus* (MRSA): burden of disease and control challenges in Europe. *Eur. Surveill*. 15:19688.
- Köck, R., Becker, K., Cookson, B., van Gemert-Pijnen, J. E., Harbarth, S., Kluytmans, J., et al. (2014). Systematic literature analysis and review of targeted preventive measures to limit healthcareassociated infections by meticillin-resistant *Staphylococcus aureus*. *Eur. Surveill.* 19:20860. doi: 10.2807/1560-7917.es2014.19.29. 20860
- Murray, J., Evans, P., and Hukins, D. W. (1980). Light-scattering methods for antibiotic sensitivity tests. J. Clin. Pathol. 33, 995–1001. doi: 10.1136/jcp.33. 10.995
- Phillips, D. T. (1971). Evolution of a light scattering photometer. *BioScience* 21, 865–867. doi: 10.1063/1.4821175
- Poulsen, M. O., Jacobsen, K., Thorsing, M., Kristensen, N. R., Clasen, J., Lillebaek, E. M., et al. (2013). Thioridazine potentiates the effect of a beta-lactam antibiotic against *Staphylococcus aureus* independently of *mecA* expression. *Res. Microbiol.* 164, 181–188. doi: 10.1016/j.resmic.2012. 10.007
- Roca, I., Akova, M., Baquero, F., Carlet, J., Cavaleri, M., Coenen, S., et al. (2015). The global threat of antimicrobial resistance: science for intervention. *New Microbes New Infect.* 6, 22–29. doi: 10.1016/j.nmni.2015. 02.007
- Sabath, L. D., Garner, C., Wilcox, C., and Finland, M. (1975). Effect of inoculum and of beta-lactamase on the anti-staphylococcal activity of thirteen penicillins and cephalosporins. *Antimicrob. Agents Chemother.* 8, 344–349. doi: 10.1128/ AAC.8.3.344

- Stull, V. R. (1973). Clinical laboratory use of differential light scattering. I. Antibiotic susceptibility testing. *Clin. Chem.* 19, 883–890.
- van Belkum, A., Durand, G., Peyret, M., Chatellier, S., Zambardi, G., Schrenzel, J., et al. (2013). Rapid clinical bacteriology and its future impact. *Ann. Lab. Med.* 33, 14–27. doi: 10.3343/alm.2013.33.1.14
- Wyatt, P. J. (1972). Light scattering in the microbial world. J. Colloid Interface Sci. 39, 479–491. doi: 10.1016/0021-9797(72)90056-2

**Conflict of Interest Statement:** EI received congress travel support by BacterioScan Inc.

The other authors declare that the research was conducted in the absence of any commercial or financial relationships that could be construed as a potential conflict of interest.

Copyright © 2017 Idelevich, Hoy, Görlich, Knaack, Grünastel, Peters, Borowski and Becker. This is an open-access article distributed under the terms of the Creative Commons Attribution License (CC BY). The use, distribution or reproduction in other forums is permitted, provided the original author(s) or licensor are credited and that the original publication in this journal is cited, in accordance with accepted academic practice. No use, distribution or reproduction is permitted which does not comply with these terms.